## DATASHEET

The Child Airway Management Trainer is perfect for practicing airway management skills like intubation, ventilation and suctioning techniques on pediatric patients.

It represents a head of a 5-year-old child featuring realistic size and anatomy including mouth, tongue, oral pharynx, larynx, epiglottis, vocal cords, trachea, and esophagus.

The separate lungs allow auscultation and the stomach bladder indicates esophageal insufflation.

The trainer allows you to practice oral, digital, and nasal intubation, as well as ET tubes, EOA, EGTA, PTL, laryngeal masks, Combitube® and King Systems insertion.

## Child Airway <u>Management Tr</u>ainer



## **KEY BENEFITS**

- Simulating a 5-year-old child head
- Attachable to the Child CRiSis™ Manikin
- Airway management training
- Jaw-thrust maneuver
- Suctioning techniques
- Carotid pulse, manual
- With lungs and stomach

## ORDERING SPECIFICATIONS

Item no.

LF03610U

**Child Airway Management Trainer** incl. head, chest plate, lungs, stomach, pulse bulb, spray lubricant



Ambu A/S Baltorpbakken 13 2750 Ballerup Denmark T +45 72 25 20 00 F +45 72 25 20 50 ambu.com