











Ambu® aScope Broncho	No	No	
D			
Videolaringoscopio King Vision	No	No	
Videolaringoscopio King™ Vision aBlade	No	No	
King Vision cable video	No	No	
VivaSight-DL	No	No	
VivaSight-SL	No	No	
VivaSight-EB	No	No	
Anaesthesia			
Dual Limb Breathing Circuits	Yes, circuit kit contains part/-s with DEHP (CASRN 117-81-7) if symbol is on label	No	
King F2	Yes, circuit kit contains part/-s with DEHP (CASRN 117-81-7) if symbol is on label	No	
King Flex2	Yes, circuit kit contains part/-s with DEHP (CASRN 117-81-7) if symbol is on label	No	
King F Flex2	No	No	
King Universal F	No	No	
King Universal F Flex	Yes, circuit kit contains part/-s with DEHP (CASRN 117-81-7) if symbol is on label	No	
King Ped F2 and King Flex 2 - Pediatric Breathing Circuits	Yes, circuit kit contains part/-s with DEHP (CASRN 117-81-7) if symbol is on label	No	
King Jackson Rees and King Mapleson	No	No	
King Transport and PACU Kits	No	No	
King Breathing Bags (STD and colored)	No	No	
HMEs and Filters	No	No	
Xenon	No	No	
Sampling Line	Yes, Contain DEHP (CASRN 117-81-7)	No	
King LT-D	Yes, Contain DEHP (CASRN 117-81-7)	No	
King LTS-D	No	No	

Producto	Contenido de SEP clasificadas por encima del 0.1% (p / p)	Látex de caucho natural	Símbolo en la etiqueta de producto
Nitric Oxide	No	No	
Vent Circuits	No	No	
Extension Sets	No	No	
Flex Connectors	No	No	
Headstraps	No	No	
Tube Tree	No	No	
Manifolds	No	No	
Ambu® Disposable Face mask	No	No	
Ambu® Open Cuff Silicone Mask	No	No	
Ambu® Silicone Face Mask	No	No	
Ambu® Transparent Face Mask	No	Sí	
Ambu® UltraSeal	No	No	
Ambu® King Mask	No	No	
Ambu® AuraOnce™	No	No	
Ambu® AuraFlex™	No	No	
Ambu® AuraStraight™	No	No	
Ambu® Aura-i™	No	No	
Ambu® Aura40™	No	No	
Ambu® AuraGain	No	No	
Ambu® Cuff Pressure Gauge	No	No	
Emergency Care			
Ambu® SPUR II	No	No	
Ambu® Mark IV	No	No	
Ambu® Mark IV Baby	No	No	
Ambu® Military Mark III	No	No	
Ambu® Oval Silicone resuscitator	No	No	
Ambu® Oval Plus Silicone resuscitator	No	No	
Ambu® PEEP Valves	No	No	
Ambu® Disposable Pressure Manometer	No	No	

Producto	Contenido de SEP clasificadas por encima del 0.1% (p / p)	Látex de caucho natural	Symbol on the product label
Ambu® TubeCheck	Sí, contain DEHP (CASRN 117-81-7)	No	
Ambu® ET Tube holder	No	No	
Ambu® LifeKey	No	Sí*	
Ambu® ResKey	No	Sí*	
Ambu® Res-Cue Mask	No	No	
Ambu® Res-Cue Pump	No	No	
Ambu® Twin Pump	No	No	
Ambu® Uni-Suction Pump	No	No	
Ambu® Perfit ACE	No	No	
Ambu® Mini Perfit ACE	No	No	
Ambu® Perfit	No	No	
Ambu® Redi-ACE Adult	No	No	
Ambu® Redi-ACE Mini	No	No	
Ambu® Head Wedge	No	No	

*No tiene contacto directo con el paciente la parte que contiene látex de caucho natural.