

Ambu® King Mask

Mascarilla Facial Desechable

HOJA TÉCNICA

The Ambu King Mask is an anatomically shaped, disposable anaesthetic facemask that is designed to support a tight seal around the facial anatomy with minimal pressure to the face. The Ambu® King Mask has a thin cushion which improves the ability to create a tight seal on a variety of patients. Furthermore, by inflating or deflating the cuff using the check valve in order to suit the patient's anatomical contours, the seal can be adjusted.

The Ambu King Mask is listed as MR Safe and is therefore suitable for use in an MR environment.

The Ambu® King Mask is also available in a scented version with Fresh Scent, Bubble Gum, Strawberry and Cherry scents. This can have a calming effect during induction, especially for children¹.

Cinco aromas

- Sin perfume (Azul)
- Menta (Verde)
- Chicle (Morado)
- Fresa (Rosa)
- Cereza (Rojo)

Siete tallas

- Talla 1 Neonato
- Talla 2 Bebés
- Talla 3 Niños de 1 a 3 años
- Talla 4 Niño/Adulto pequeño
- Talla 5 Adulto mediano
- Talla 6 Adulto grande
- Talla 7 RD Adulto extra grande



CARACTERÍSTICAS PRINCIPALES

Seal can be enhanced by deflating or inflating the cushion via the check valve

Cushion is thin and pliable, for an excellent face seal

Ribbed on the edge of the cushion, for improved grip

Anatomically shaped cuff for improved adjustment to the face with minimal pressure

Clear dome for easy observation of the patient's condition

MR Safe

Available with or without hook ring*

¹ Fukumoto M, Arima H, Ito S, Takeuchi N, Nakano H. Distorted perception of smell by volatile agents facilitated inhalational induction of anesthesia. Paediatr Anaesth. 2005 Feb; 15 (2): 98-101.

*No disponible para tamaños 1 y 2.

Ambu

