

DATASHEET

The Sani-Baby CPR Manikin is an economical infant manikin and features realistic anatomical landmarks including sternum, rib cage and substernal notch.

The manikin allows the training of the BLS algorithm as well as the realistic head tilt to open the airways. A successful ventilation will be indicated by a realistic chest rise.

Sanitary System

For preventing cross-contamination, a face shield/lung system is inserted easy and quickly into the manikin. This system allows a realistic expiration and no cleaning or disinfecting is necessary.

Sani-Baby CPR



KEY BENEFITS

- Economical infant manikin for BLS training
- Realistic face and molded hair
- Anatomically correct with palpable landmarks
- Visible chest rise during ventilation
- Sanitary system prevents cross-contamination: single-use face shield/lung system

ORDERING SPECIFICATIONS

Item no.

100-2121U	Sani-Baby CPR Manikin Full-body manikin with diaper, 3 face shield/lung systems
100-2125U	Replacement Face Shield/Lung System for Sani-Baby (set of 100)

Ambu

📍 Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com